



Product Datasheet

rFel d 1 (Felis domesticus 1.0101) (orb1088525)


Description	Recombinant protein rFel d 1 is expressed in Escherichia coli. DNA sequence encoding 175 AAs was fus...
Tested Applications	ELISA, FC, WB
Storage	Store at -20°C to -80°C. Reconstitute in sterile deionized water. Use reconstituted product immediately or aliquot for further storage at -20°C to -80°C.
Note	For research use only
Purity	Purified by sequential steps of ion exchange and affinity chromatography. Lyophilized from a 0.2 µm filtered solution in Storage buffer (100 mM Tris, 150 mM NaCl, 0.01% Tween 20, pH 8.0).
Uniprot ID	P30438, P30440
Expiration Date	12 months from date of receipt.



Recombinant allergen rFel d 1 purity ver...



ELISA (enzyme-linked immunosorbent assay...



Flow cytometry dot-plot staining pattern...