

## Product Datasheet

# Recombinant Human VEGFA Protein, C-His (orb1088542)

<b>Description</b>	Recombinant Human Vascular Endothelial Growth Factor (VEGF165) protein is produced by Mammalian expr...
<b>Tested Applications</b>	CLIA, EIA, LF
<b>Target</b>	VEGFA
<b>Preservatives</b>	PBS,20%Glycerol,pH 7.4
<b>Form/Appearance</b>	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4
<b>Concentration</b>	Lot specific(±20%)
<b>Storage</b>	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
<b>Tag</b>	His-tag
<b>Note</b>	For research use only
<b>Application notes</b>	Use: Calibration/quality control. Quality control: Mol Mass: 19.8 kDa AP Mol Mass: 24-27 kDa Purity: ≥90% as determined by reducing SDS-PAGE
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Source</b>	HEK293
<b>MW</b>	20.8 kDa
<b>Uniprot ID</b>	<b>P15692</b>
<b>Expiration Date</b>	6 months from date of receipt.