

Product Datasheet

Recombinant Human ADP Protein (orb1088544)

Description	Recombinant Human adiponectin (ADP) protein is produced by Mammalian expression system and the target...
Tested Applications	CLIA, EIA, LF
Target	ADP
Preservatives	PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4
Concentration	Lot specific(±20%)
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Tag	His-tag
Note	For research use only
Application notes	Use: Calibration/quality control. Quality control: Mol Mass: 26.4 kDa AP Mol Mass: 28.5 kDa Purity: ≥90% as determined by reducing SDS-PAGE
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293
MW	30 kDa
Uniprot ID	Q15848
Expiration Date	6 months from date of receipt.