

Product Datasheet

Recombinant Human Lp-PLA2 (C-His) (orb1088563)

Description	Recombinant Human Lipoprotein-associated phospholipase A2(LP-PLA2) protein is produced by E.coli exp...
Tested Applications	CLIA, EIA, LF
Target	Lp-PLA2
Preservatives	20 mM Tris, 150 mM NaCl, 10% Glycerol, pH 9.0.
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20 mM PB, 150 mM NaCl, 5% glycerol, pH 7.4
Concentration	Lot specific(±20%)
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Tag	His-tag
Note	For research use only
Application notes	Use: Calibration/quality control. Quality control: Mol Mass: 47.85 kDa AP Mol Mass: 47-50 kDa Purity: ≥90% as determined by reducing SDS-PAGE
Purity	> 85% as determined by SDS-PAGE.
Source	E.coli
MW	48.7 kDa
Expiration Date	6 months from date of receipt.