

Product Datasheet

Recombinant Human IL6 (C-His) (orb1088571)

Description	Recombinant Human Interleukin 6 protein is produced by Mammalian expression system and the target ge...
Tested Applications	CLIA, EIA, LF
Target	IL-6
Preservatives	20 mM PB, 150 mM NaCl, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4
Concentration	Lot specific(±20%)
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Tag	His-tag
Note	For research use only
Application notes	Use: Calibration/quality control. Quality control: Mol Mass: 21.67 kDa AP Mol Mass: 24 kDa Purity: ≥90% as determined by reducing SDS-PAGE
Purity	> 95% as determined by SDS-PAGE
Source	HEK293
Expiration Date	6 months from date of receipt.