

Product Datasheet

phospho-RACGAP1 (Ser387) Polyclonal Antibody (orb1088659)

Description	phospho-RACGAP1 (Ser387) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human RACGAP1 around the phosphorylation site of Ser387:RI(p-S)GC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	71kDa
Uniprot ID	Q9H0H5
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.