

Product Datasheet

Retinoic acid induced 17 Polyclonal Antibody (orb1088665)

Description	Retinoic acid induced 17 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Gallus, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Retinoic acid induced 17:51-150/1067
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	115kDa
Uniprot ID	Q9ULJ6
Dilution Range	IHC-P=1:100-500 IHC-F=1:100-500 Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.