


Product Datasheet

Protein Phosphatase 5 antibody (orb108868)

Description	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	This antibody can be used for Western blotting (5-10 µg/ml). Other applications being evaluated. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Source	Rabbits were immunized with a Synthetic peptide derived from the human protein phosphatase 5 proteins
Uniprot ID	P53041
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 and MCF-7 ...