

Product Datasheet

phospho-TNIK (Ser769) Polyclonal Antibody (orb1088734)

Description	phospho-TNIK (Ser769) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TNIK around the phosphorylation site of Ser769:EG(p-S)PV
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	150kDa
Uniprot ID	Q9UKE5
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.