

## Product Datasheet

# phospho-HMGCR (Ser872) Polyclonal Antibody (orb1088770)

<b>Description</b>	phospho-HMGCR (Ser872) Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Canine, Gallus, Human, Mouse, Porcine, Rabbit, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human HMGCR around the phosphorylation site of Ser872:NR(p-S)KI
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -55°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~97 kDa
<b>Uniprot ID</b>	<b>P04035</b>
<b>Dilution Range</b>	Flow-Cyt=1ug/Test
<b>Expiration Date</b>	12 months from date of receipt.