

Product Datasheet

PAFR/PAF Receptor Polyclonal Antibody (orb1088872)

Description	PAFR/PAF Receptor Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Gallus, Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PAFR:231-342/342
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	38kD
Uniprot ID	P25105
Dilution Range	WB=1:500-2000 IHC-P=1:100-500
Expiration Date	12 months from date of receipt.