


Product Datasheet

DAP12 antibody (orb108888)

Description	Rabbit polyclonal to DAP12
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for immunoprecipitation (2-4 µg/mg of protein lysate). Optimal concentration should be evaluated by serial dilutions.
Source	Rabbit immunized with recombinant DAP12 protein
Uniprot ID	O43914
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of NKp44 po...