

Product Datasheet

FEZ2 Polyclonal Antibody (orb1088945)

Description	FEZ2 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Equine, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human FEZ2/Zygin2:165-240/353
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	46 kD
Uniprot ID	Q9UHY8
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.