

Product Datasheet

Phospho-Progesterone Receptor(Ser190) Polyclonal Antibody (orb1089023)

Description	Phospho-Progesterone Receptor(Ser190) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat Progesterone Receptor around the phosphorylation site of Ser190
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	110 kD/ 120kD
Uniprot ID	Q63449
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.