

Product Datasheet

Acetyl Coenzyme A Carboxylase Polyclonal Antibody (orb1089038)

Description	Acetyl Coenzyme A Carboxylase Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Acetyl Coenzyme A Carboxylase 1/2:2307-2355/2458
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	265kD
Uniprot ID	O00763
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.