

## Product Datasheet

### Fbx32 Polyclonal Antibody (orb1089047)

<b>Description</b>	Fbx32 Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human MAFbx:31-130/355
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	42kD
<b>Uniprot ID</b>	<a href="#">Q969P5</a>
<b>Dilution Range</b>	Flow-Cyt=2ug/Test
<b>Expiration Date</b>	12 months from date of receipt.