

## Product Datasheet

### Nano-Tag (9) monoclonal antibody (orb1089152)

<b>Description</b>	Nano-Tag (9) monoclonal antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH conjugated synthetic peptide Nano-Tag (9)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB=1:1000-10000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Expiration Date</b>	12 months from date of receipt.