

Product Datasheet

Avi-Tag monoclonal antibody (orb1089162)

| | |
|----------------------------|---|
| Description | Avi-Tag monoclonal antibody |
| Species/Host | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Recombinant Avi-Tag(GLNDIFEAQKIEWHE) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Dilution Range | WB=1:1000-5000 |
| Expiration Date | 12 months from date of receipt. |