

Product Datasheet

HSV-Tag monoclonal antibody (orb1089165)

Description	HSV-Tag monoclonal antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide HSV-Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB=1:1000-5000
Expiration Date	12 months from date of receipt.