

Product Datasheet

RGC32/C13orf15 Polyclonal Antibody (orb1089207)

Description	RGC32/C13orf15 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Gallus, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human RGC32:65-137/137
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	15kDa
Uniprot ID	Q9H4X1
Dilution Range	Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.