

Product Datasheet

phospho-CK II beta (Ser209) Polyclonal Antibody (orb1089245)

Description	phospho-CK II beta (Ser209) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CK II beta around the phosphorylation site of Ser209:FK(p-S)PV
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	24kDa
Uniprot ID	P67871
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.