

Product Datasheet

Adiponectin Receptor 2 antibody (orb108925)

Description	Rabbit polyclonal antibody to Adiponectin Receptor 2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for Western blotting (5-10µg/ml) and ELISA. Other applications not yet tested. Cross reacts with Human ADR1. Optimal concentration should be evaluated by serial dilutions.
Source	Rabbits were immunized with a synthetic peptide near the C terminus of human adiponectin receptor 2 proteins
Uniprot ID	Q86V24
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.