

Product Datasheet

Akt (Ab-450) polyclonal antibody [KU] (orb1089265)

| | |
|----------------------------|---|
| Description | Akt (Ab-450) polyclonal antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Immunogen | Synthetic peptide, corresponding to Human Akt. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 56 kDa |
| Uniprot ID | Q9Y243 , P31749 , P31751 |
| Dilution Range | WB: 1:1000~1:2000 IHC: 1:50~1:200 IF: 1:50~1:200 |
| Expiration Date | 12 months from date of receipt. |