

Product Datasheet

Akt(Ab-450) polyclonal antibody \square KU \square (orb1089273)

Description	Akt(Ab-450) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide, corresponding to Human Akt.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 60 kDa
Uniprot ID	Q9Y243 , P31749 , P31751
Dilution Range	WB: 1:2000~1:5000 IF: 1:50~1:200
Expiration Date	12 months from date of receipt.