

## Product Datasheet

### Centrinone (orb1089329)

<b>Description</b>	Centrinone is a selective and reversible inhibitor of polo-like kinase 4 (Plk4), a serine-threonine ...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{26}H_{25}F_2N_7O_6S_2$
<b>MW</b>	633.65
<b>CAS Number</b>	[1798871-30-3]
<b>SMILES</b>	<chem>CC1=CC(=NN1)NC2=NC(=NC(=C2OC)N3CCOCC3)SC4=C(C=C(C=C4)S(=O)(=O)CC5=C(C=CC=C5)[N+](=O)[O-])F)F</chem>
<b>Expiration Date</b>	12 months from date of receipt.