

## Product Datasheet

### NSC45573 (orb1089366)

<b>Description</b>	Disperse Fast Yellow 4K is a color additive used in dyes.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{19}H_{16}N_4O$
<b>MW</b>	316.35654
<b>CAS Number</b>	[6300-37-4]
<b>SMILES</b>	<chem>OC1=C(C)C=C(C=C1)\N=N\C2=CC=C(\N=N\C3=CC=CC=C3)C=C2</chem>
<b>Expiration Date</b>	12 months from date of receipt.