

Product Datasheet

LY-3475070 (orb1089369)

Description	LY-3475070 is a potent and selective CD73 Inhibitor with IC50 of 28 nM. LY3475070 Alone or in Combin...
Note	For research use only
Formula	$C_{15}H_{18}N_4O_2$
MW	286.335
CAS Number	[2375815-63-5]
SMILES	<chem>O=C1NC(C(C2=NN=C(C)C([C@@H]3[C@@H](C(C)C)C3)=C2)=CN1)=O</chem>
Expiration Date	12 months from date of receipt.