

Product Datasheet

3M-011 (orb1089456)

Description	3M-011 is a potent TLR7/8 agonist and cytokine inducer; stimulates type I interferon (IFN) and other...
Note	For research use only
Formula	$C_{18}H_{25}N_5O_3S$
MW	391.49
CAS Number	[642473-62-9]
SMILES	<chem>CS(=O)(NC(C)(C)CN1C(COCC)=NC2=C1C3=CC=CC=C3N=C2N)=O</chem>
Expiration Date	12 months from date of receipt.