

## Product Datasheet

### KUN56321 (orb1089486)

<b>Description</b>	KUN56321 is a luminescent agent. The application of KUN56321 comprises: coating the film of KUN56321...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{23}H_{15}BRN_2$
<b>MW</b>	399.28
<b>CAS Number</b>	[1771756-32-1]
<b>SMILES</b>	<chem>C1C=C(C=CC=1)N1C(=NC2C3C(=CC=CC=3)C=CC1=2)C1C=CC(=CC=1)Br</chem>
<b>Expiration Date</b>	12 months from date of receipt.