

Product Datasheet

Val-Lys(Boc)-PAB (orb1089496)

Description	Val-Lys(Boc)-PAB is a ADC linker. Val-Lys(Boc)-PAB was used to prepare camptothecin peptide conjugat...
Note	For research use only
Formula	$C_{23}H_{38}N_4O_5$
MW	450.58
CAS Number	[1432969-86-2]
SMILES	<chem>CC(C)(OC(NCCCC[C@@H](C(NC1=CC=C(CO)C=C1)=O)NC([C@H](C(C)C)N)=O)=O)C</chem>
Expiration Date	12 months from date of receipt.