

## Product Datasheet

### PI3Kδ-IN-1 (orb1089497)

<b>Description</b>	PI3Kδ-IN-1 is a potent, selective, and efficacious PI3Kδ inhibitor with an IC50 of 1.7 nM.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>22</sub> H <sub>20</sub> F <sub>3</sub> N <sub>7</sub> O <sub>2</sub>
<b>MW</b>	471.4351
<b>CAS Number</b>	[1911564-39-0]
<b>SMILES</b>	<chem>FC1C([H])=C(C2C([H])=C([H])C(C#N)=C(C=2[H])N2C(C(C([H])([H])[H])(C([H])([H])[H])N(C(C([H])([H])[H])=O)C([H])([H])C2([H])[H])=O)N2C=1C(N([H])[H])=NC([H])=N2)(F)F</chem>
<b>Expiration Date</b>	12 months from date of receipt.