

---

## Product Datasheet

### IL13 antibody (orb10895)

**Description**

IL13 antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

**Applications**
**Immunogen**

KLH conjugated synthetic peptide derived from human IL13. Please contact us for the exact immunogen sequence. The peptide is available as orb374729.

**Target**

IL13

**Form/Appearance**

10 mM PBS, 0.02% sodium azide

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**MW**


12 kDa

**Uniprot ID**
[P20109](#), [P35225](#), [P42203](#)
**NCBI**
[002188](#), [21](#), [002179](#)
**Entrez**
[3596](#)
**Dilution Range**


IHC-P: 1:50-300 (based on 0.5 mg/ml)

**Expiration Date**


12 months from date of receipt.



IHC-P image of rat epididymis tissue using...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of paraffin...