

## Product Datasheet

### Donafenib (Sorafenib D3) (orb1089591)

<b>Description</b>	Sorafenib D3 (Bay 43-9006 D3) is the deuterium labeled Sorafenib. Sorafenib is a multikinase inhibit...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{16}ClF_3N_4O_3$
<b>MW</b>	467.8434391
<b>CAS Number</b>	[1130115-44-4]
<b>SMILES</b>	<chem>[2H]C(NC(C1=NC=CC(OC2C=CC(NC(NC3C=CC(Cl)=C(C(F)(F)F)C=3)=O)=CC=2)=C1)=O)([2H])[2H]</chem>
<b>Expiration Date</b>	12 months from date of receipt.