

Product Datasheet

Heparin sodium salt (MW 15kDa) (orb1089611)

Description	Heparin sodium salt (MW 15kDa) (Sodium heparin (MW 15kDa)) is a polymer of Heparin with the molecula...
Note	For research use only
Formula	$C_{12}H_{17}NO_{20}S_3$
MW	591.453081846237
CAS Number	[9041-08-1]
SMILES	<chem>O[C@H]1[C@H](O[C@H]2[C@H](OS(O)(=O)=O)[C@@H](O)[C@H](OC)[C@H](C(O)=O)O2)[C@@H](COS(O)(=O)=O)O[C@H](OC)[C@@H]1NS(O)(=O)=O.[n].[xNa]</chem>
Expiration Date	12 months from date of receipt.