

## Product Datasheet

### CAN-508 (orb1089674)

<b>Description</b>	CAN-508 is a potent, selective CDK9 inhibitor with IC50 of 0.35 uM (CDK9/cyclin T1), displays > 35 f...
<b>Note</b>	For research use only
<b>Formula</b>	$C_9H_{10}N_6O$
<b>MW</b>	218.22
<b>CAS Number</b>	[140651-18-9]
<b>SMILES</b>	<chem>OC1=CC=C(C=C1)\N=N\C2C(N)=NNC=2N</chem>
<b>Expiration Date</b>	12 months from date of receipt.