

## Product Datasheet

### KT-109 (orb1089682)

<b>Description</b>	KT-109 is a potent, selective inhibitor of DAGLβ with IC50 of 42 nM, displays about 60-fold selecti...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{27}H_{26}N_4O$
<b>MW</b>	422.5215459
<b>CAS Number</b>	[1402612-55-8]
<b>SMILES</b>	<chem>C1(C2C=CC=CC=2)C=CC(C2N=NN(C(N3CCCCC3CC3C=CC=CC=3)=O)C=2)=CC=1</chem>
<b>Expiration Date</b>	12 months from date of receipt.