

## Product Datasheet

### LY2033298 (orb1089707)

<b>Description</b>	LY2033298 is a robust allosteric potentiator that is highly selective for the human M4 muscarinic ac...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{13}H_{14}ClN_3O_2S$
<b>MW</b>	311.7871604
<b>CAS Number</b>	[886047-13-8]
<b>SMILES</b>	<chem>NC1C2C(=NC(OC)=C(CI)C=2C)SC=1C(NC1CC1)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.