

Product Datasheet

Propargyl-PEG2-acid (orb1089767)

Description	Propargyl-PEG2-acid is a non-cleavable 2 unit PEG ADC linker used in the synthesis of antibody-drug ...
Note	For research use only
Formula	$C_8O_4H_{12}$
MW	172.1785
CAS Number	[1859379-85-3]
SMILES	<chem>O(C([H])([H])C([H])([H])C(=O)O[H])C([H])([H])C([H])([H])OC([H])([H])C#C[H]</chem>
Expiration Date	12 months from date of receipt.