

## Product Datasheet

### L-798106 (orb1089778)

<b>Description</b>	L-798106 is potent and highly selective prostanoid EP3 receptor antagonist (K <sub>i</sub> =0.3 nM), it also has ...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>27</sub> H <sub>22</sub> BRNO <sub>4</sub> S
<b>MW</b>	536.4368854
<b>CAS Number</b>	[244101-02-8]
<b>SMILES</b>	<chem>COC1C=CC(Br)=CC=1S(NC(/C=C\C1=CC=CC=C1CC1C=CC2=CC=CC=C2C=1)=O)(=O)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.