

## Product Datasheet

### m-PEG3-Amine (orb1089793)

<b>Description</b>	m-PEG3-Amine is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs).
<b>Note</b>	For research use only
<b>Formula</b>	$C_7H_{17}NO_3$
<b>MW</b>	163.2147825
<b>CAS Number</b>	[74654-07-2]
<b>SMILES</b>	COCCOCCOCCN
<b>Expiration Date</b>	12 months from date of receipt.