

## Product Datasheet

### D609 (orb1089804)

<b>Description</b>	D609 is a selective competitive phosphatidyl choline-specific phospholipase C (PC-PLC) inhibitor (Ki...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{11}H_{16}O_5 \cdot K$
<b>MW</b>	267.47244
<b>CAS Number</b>	[83373-60-8]
<b>SMILES</b>	[K].SC(=S)OC1C2C3C(CCC3)C(C1)C2
<b>Expiration Date</b>	12 months from date of receipt.