

Product Datasheet

Metarrestin (orb1089807)

Description	Metarrestin (ML246) is a specific inhibitor of perinuclear compartment (PNC), disrupts PNCs in PC3M-...
Note	For research use only
Formula	$C_{31}H_{30}N_4O$
MW	474.596107
CAS Number	[1443414-10-5]
SMILES	<chem>O([H])C1([H])C([H])([H])C([H])([H])C([H])(C([H])([H])C1([H])([H])N1C([H])=NC2=C(C(C3C([H])=C([H])C([H])=C([H])C=3[H])=C(C3C([H])=C([H])C([H])=C([H])C=3[H])N2C([H])([H])C2C([H])=C([H])C([H])=C([H])C=2[H])/C/1=N/[H]</chem>
Expiration Date	12 months from date of receipt.