

Product Datasheet

TM5007 (orb1089811)

Description	TM5007 is a potent inhibitor plasminogen activator inhibitor-1 (PAI-1) with IC50 of 29 uM, exhibits n...
Note	For research use only
Formula	$C_{24}H_{20}N_2O_6S_4$
MW	560.69
CAS Number	[342595-05-5]
SMILES	<chem>C(NC1SC=C(C2SC=CC=2)C=1C(O)=O)(=O)CCCCC(NC1SC=C(C2SC=CC=2)C=1C(O)=O)=O</chem>
Expiration Date	12 months from date of receipt.