

## Product Datasheet

### ML-299 (orb1089813)

<b>Description</b>	ML299 is a potent, CNS penetrant, dual phospholipase PLD1/PLD2 probe with IC50 of 6 nM/20 nM, respec...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{23}H_{26}BrFN_4O_2$
<b>MW</b>	489.389
<b>CAS Number</b>	[1426916-00-8]
<b>SMILES</b>	<chem>BrC1C=CC(=CC=1)C(N[C@@H](C)CN1CCC2(C(NCN2C2C=CC=C(C=2)F)=O)CC1)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.