

## Product Datasheet

### XS018661 (orb1089839)

<b>Description</b>	The first-in-class dual inhibitor of ENL (Kd = 754 ± 45 nM) and its paralog AF9 (Kd = 523 ± 53 nM)
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>22</sub> H <sub>24</sub> N <sub>4</sub> O <sub>3</sub>
<b>MW</b>	392.46
<b>SMILES</b>	<chem>O=C(OC)C1=CC=C(C(C(NC2=CC=C3NC(CN4CCCCC4)=NC3=C2)=O)C=C</chem> 1
<b>Expiration Date</b>	12 months from date of receipt.