

Product Datasheet

G244-LM (orb1089903)

Description	G244-LM is a potent and specific inhibitor of tankyrase 1/2 that inhibits Wnt signaling.
Note	For research use only
Formula	$C_{18}H_{22}N_4O_3S_2$
MW	406.52
CAS Number	[1563007-08-8]
SMILES	<chem>C1(N2CCN(C3=CC=CC=C3S(C)(=O)=O)CC2)NC(=O)C2CSCCC=2N=1</chem>
Expiration Date	12 months from date of receipt.