

Product Datasheet

FOXP2 antibody (orb1089919)

| | |
|----------------------------|---|
| Description | Human polyclonal antibody to FOXP2 |
| Species/Host | Human |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, IP |
| Immunogen | This antibody was obtained by recombinant antibody (rAb) phage display recognizing FOXP2 protein under non-denaturing conditions. |
| Target | FOXP2 |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | 1 mg/ml |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Application notes | This antibody has been used and recommended for the following methods: ELISA; spiked-IP; FC. |
| Isotype | Fab fragment |
| Clonality | Monoclonal |
| Purity | Purified |
| Uniprot ID | O15409 |
| Expiration Date | 12 months from date of receipt. |