

Product Datasheet

PRDM9 antibody (orb1089927)

Description	Human polyclonal antibody to PRDM9
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing PRDM9 protein under non-denaturing conditions.
Target	PRDM9
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has been used and recommended for the following methods: ELISA; IF; IP-MS; ChIP.
Isotype	Fab fragment
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q9NQV7
Expiration Date	12 months from date of receipt.