

Product Datasheet

SUV420H2 antibody (orb1089945)

Description	Rabbit polyclonal antibody to SUV420H2
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing SUV420H2 protein under non-denaturing conditions.
Target	SUV420H2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has been used and recommended for the following methods: ELISA; IP-MS; ChIP.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q86Y97
Expiration Date	12 months from date of receipt.