

## Product Datasheet

### L3MBTL1 antibody (orb1089952)

<b>Description</b>	Rabbit polyclonal antibody to L3MBTL1
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ChIP, ELISA, FC, IF, IP
<b>Immunogen</b>	This antibody was obtained by recombinant antibody (rAb) phage display recognizing L3MBTL1 protein under non-denaturing conditions.
<b>Target</b>	L3MBTL1
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody has been used and recommended for the following methods: ELISA; IF; spiked IP; IP-MS; ChIP; FC.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">Q9Y468</a>
<b>Expiration Date</b>	12 months from date of receipt.